



## Anti-CEP-1 ELISA (IgG)



- High specificity for RA, thus a supplementary serological parameter for securing a diagnosis
- Research into pathogenesis of RA and subtyping of the disease
- Potentially relevant for monitoring treatment responses and assessing therapy success

### Technical data

<b>Antigen</b>	Citrullinated $\alpha$ -enolase peptide 1 (CEP-1)
<b>Calibration</b>	Quantitative, in relative units per milliliter (RU/ml)
	Calibrator 1: 200 RU/ml
	Calibrator 2: 20 RU/ml
	Calibrator 3: 2 RU/ml
	Recommended upper threshold value for normal range (cut-off): 20 RU/ml
<b>Sample dilution</b>	Serum or plasma; 1:101 in sample buffer
<b>Reagents</b>	Ready for use, with the exception of the wash buffer (10x)
<b>Test procedure</b>	30 min / 30 min / 15 min. Room temperature. Fully automatable
<b>Measurement</b>	450 nm, Reference wavelength between 620 nm and 650 nm
<b>Kit format</b>	96 break-off wells. Kit includes all necessary reagents.
<b>Order number</b>	EA 151b-9601 G

### Clinical significance

In suspected cases of rheumatoid arthritis (RA), mainly rheumatoid factors (RF) and autoantibodies against cyclic citrullinated peptides (CCP) are serologically determined alongside general inflammation parameters. RF are found in 70 to 90% of sera from RA patients, but they also occur in patients with other autoimmune and infectious diseases. The detection of autoantibodies against CCP which is optimised for high specificity (97%) and sensitivity (72%) is on the other hand an established standard test. However, although available tests for the detection of autoantibodies against citrullinated protein/peptide antigens (ACPA) such as CCP allow the detection of ACPA, they do not deliver any information about the nature of the autoantibodies. To gain more information, the reactivity of the citrullinated proteins that actually occur in the inflamed joints of RA patients, for example fibrinogen/fibrin, vimentin, collagen type II, fibronectin and  $\alpha$ -enolase, must be investigated. Citrullinated  $\alpha$ -enolase was discovered as an autoantigen in RA in 2005. The epitope within the  $\alpha$ -enolase "citrullinated  $\alpha$ -enolase peptide 1" (CEP-1) which is responsible for the autoimmune reaction was described as the relevant autoantigen in 37 to 62% of RA patients by Prof. Lundberg. Antibodies against this peptide show a high specificity, and are only found in the sera of 2 to 3% of healthy blood donors or control patients.

### Diagnostic application

Antibodies against CEP-1 are a supplementary parameter for diagnosis of RA and for research into the pathogenesis and subtyping of the disease. In collaboration with Prof. Venables (Kennedy Institute of Rheumatology, University of Oxford), a commercial CE-certified Anti-CEP-1 ELISA has been developed, for which EUROIMMUN has an exclusive license. Results from studies demonstrate a sensitivity of 50% on average and a high specificity of 98% for the test system. In 2009 it was shown in *Nature Genetics* that autoantibodies against CEP-1 are associated with a subtype of RA for which smoking is a main risk factor. Further, an association between anti-CEP-1 and *Porphyromonas gingivalis* infection, one of the main causes of periodontitis, has been reported. Since anti-CEP-1 represents an immune response to a relevant target antigen, the detection of these autoantibodies is better suited for gaining insight into the cause and pathogenesis of RA than CCP2-based test systems. Current studies are investigating whether the determination of autoantibodies against CEP-1 might also be suitable for prognosis about the response to treatment and for monitoring the disease and therapy course.



## Linearity

The linearity of the Anti-CEP-1 ELISA (IgG) was determined by assaying at least 4 serial dilutions of different patient samples. The coefficient of determination  $R^2$  was  $>0.95$  for all sera. The Anti-CEP-1 ELISA (IgG) is linear at least in the tested concentration range of 3 RU/ml to 200 RU/ml.

## Reference range

The levels of anti-CEP-1 antibodies were determined in 300 samples from apparently healthy blood donors using the EUROIMMUN ELISA. The mean concentration of antibodies against CEP-1 was 4 RU/ml with values ranging from 0.9 to 53.6 RU/ml. With a cut-off of 20 RU/ml, 1.7% of the blood donors were anti-CEP-1 positive.

Blood donors (n = 300)			
Percentile	97 <sup>th</sup>	98 <sup>th</sup>	99 <sup>th</sup>
Cut-off	16.1 RU/ml	19.5 RU/ml	23.8 RU/ml

## Reproducibility

To evaluate the test reproducibility the intra- and inter-assay coefficients of variation were determined using 3 sera. The intra-assay CVs are based on 20 measurements and the inter-assay CVs on 3 determinations repeated in 10 different test runs.

Serum	Intra-assay variation, n=20		Inter-assay variation, n=3 x 10	
	Mean value (RU/ml)	CV (%)	Mean value (RU/ml)	CV (%)
1	16.3	3.6	16.5	7.6
2	67.3	2.4	67.5	5.7
3	170.3	1.5	170.9	2.5

## Sensitivity and specificity

110 patients with rheumatoid arthritis (RA), a control panel of 238 patients with other diseases (SLE, SS, SCL, FM, SpA, PsA or infectious diseases) and 500 apparently healthy blood donors were investigated with the EUROIMMUN Anti-CEP-1 ELISA (IgG). The sensitivity of the ELISA for RA was 43.6% at a specificity of 97.6%.

Panel (n = 848)	n	Anti-CEP-1 ELISA (IgG) positive
Rheumatoid arthritis (RA)	110	48
<b>Sensitivity for RA</b>	<b>110</b>	<b>48 (43.6%)</b>
Systemic lupus erythematosus (SLE)	120	5
Sjögren's syndrome (SS)	30	1
Scleroderma (SCL)	30	2
Fibromyalgia (FM)	20	0
Axial spondyloarthritis (SpA)	20	1
Psoriatic arthritis (PsA)	2	0
Infectious diseases	16	1
Healthy blood donors	500	8
<b>Specificity for RA</b>	<b>738</b>	<b>18 (97.6%)</b>

## Literature

- Fisher BA, Plant D, Brode M, van Vollenhoven RF, Mathsson L, Symmons D, Lundberg K, Rönnelid J, Venables PJ. Antibodies to citrullinated  $\alpha$ -enolase peptide 1 and clinical and radiological outcomes in rheumatoid arthritis. *Ann Rheum Dis*. 2011 Jun;70(6):1095-8.
- EUROIMMUN AG. Stöcker W, Schlumberger W, Krüger C. In: Gressner A, Arndt T (Eds.) *Lexicon of Medical Laboratory Diagnostics*. 2. ed. Springer Medizin Verlag, Heidelberg (2012) (In German).
- Kinloch A, Tatzler V, Wait R, Peston D, Lundberg K, Donatien P, Moyes D, Taylor PC, Venables PJ. Identification of citrullinated alpha-enolase as a candidate autoantigen in rheumatoid arthritis. *Arthritis Res Ther*. 2005;7(6):R1421-9. Epub 2005 Oct 19.
- Kinloch A, Lundberg K, Wait R, Wegner N, Lim NH, Zendman AJ, Saxne T, Malmström V, Venables PJ. Synovial fluid is a site of citrullination of autoantigens in inflammatory arthritis. *Arthritis Rheum*. 2008 Aug;58(8):2287-95.
- Lundberg K, Kinloch A, Fisher BA, Wegner N, Wait R, Charles P, Mikuls TR, Venables PJ. Antibodies to citrullinated alpha-enolase peptide 1 are specific for rheumatoid arthritis and cross-react with bacterial enolase. *Arthritis Rheum*. 2008 Oct;58(10):3009-19.
- Lundberg K, Wegner N, Yucel-Lindberg T, Venables PJ. Periodontitis in RA-the citrullinated enolase connection. *Nat Rev Rheumatol*. 2010 Dec;6(12):727-30.
- Mahdi H, Fisher BA, Källberg H, Plant D, Malmström V, Rönnelid J, Charles P, Ding B, Alfredsson L, Padyukov L, Symmons DPM, Venables PJ, Klareskog L, Lundberg K. Specific interaction between genotype, smoking and autoimmunity to citrullinated  $\alpha$ -enolase in the etiology of rheumatoid arthritis. *Nature Genetics* 41 (2009) 1319-1324.